

## PowerPlay® Cold & Compression Knee Wrap


## **INSTRUCTIONS FOR USE**


2


PRODUCT #: PPKN-58


- Remove optional gel pack and place in freezer for a minimum of 1 hour.
  - NOTE: Each time the gel wrap thaws, refreeze for 1 hour.
- Place a protective barrier between the gel pack and the skin. A
  cloth sleeve is included in your package. Any thin cloth barrier
  will be effective. Remove gel from freezer, velcro the gel wrap
  on the inside of the brace.
- Place the brace over the knee and secure the two straps around the knee to achieve a snug fit.

- 4. Connect the hose to a port on the PowerPlay pump. Press the Power button to turn on the unit.
- Press the Set Pressure button associated with that port until the desired pressure appears on the display.
  - NOTE: Pressure cannot be altered during an inflation cycle.









For questions about this product, call 855-PWR-PLAY (855-797-7529)





## PowerPlay® Cold & Compression Elbow Wrap

## **INSTRUCTIONS FOR USE**

PRODUCT #: PPEL-22

- Remove optional gel pack and place in freezer for a minimum of 1 hour.
  - NOTE: Each time the gel wrap thaws, refreeze for 1 hour.
- Remove gel from freezer and secure the gel pack to the inside of the brace.
- Place a protective barrier between the gel pack and the skin. A
  cloth sleeve is included in your package. Any thin cloth barrier
  will be effective. Place the elbow onto the brace and secure
  the two straps around the arm to achieve a snug fit.
- Connect the hose to a port on the PowerPlay pump. Press the Power button to turn on the unit.
- 5. Press the Pressure Set button associated with that port until the desired pressure appears on the display.
  - NOTE: Pressure cannot be altered during an inflation cycle.

